Simport® Scientific Histology Cassettes - Instructions for use

EN

Tissue and biopsy cassettes are designed for processing and embedding of specimens for histological analysis. Cassettes are used for fixation, processing, infiltration, embedding and sectioning of fixed or fresh tissue samples. They provide a writing surface for identification and for archiving.

INSTRUCTIONS FOR USE

1. Cassette models should be specific to each tissue or biopsy size; choose the right model for each specimen size
 - Smaller specimens can be contained with the use of biopsy foam pads, biopsy bags or wrapping if necessary
 - Larger specimens may be trimmed to ensure a good fit for optimal processing. Do not overfill the cassette
2. Identify each cassette with the appropriate specimen/patient information during sampling; use identification methods that are resistant to histological solvents.
3. A separate or attached lid tightly closed will securely hold the specimen. If metal lids are used, do not use microwave techniques.
4. Process samples according to specific laboratory techniques. Leave enough space between cassettes to allow adequate fluid transfer
 - Leaving cassettes in acid or histological solvents for long periods may adversely affect design, performance and color. Metal lids can be corroded by acid, which could compromise integrity of the closure
 - Metal lids should be cleaned before reuse
 - Plastic lids are single use only.
5. When embedding tissue or biopsy samples, fill the back of the cassette thoroughly with paraffin to ensure proper support of the specimen in the microtome.
6. Avoid air pockets in hot paraffin when embedding before it solidifies.

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