EN

HISTOSETTE; MACROSETTE; MICROMESH; MICROSETTE; SLIMSETTE; SWINGSETTE; UNISETTE, SUPERFLO

Are single use. Tissue and biopsy cassettes are designed for processing and embedding, fixation, infiltration and manual sectioning of fixed or fresh tissue samples for histological analysis. They provide a writing surface for identification and for archiving as well as a support for microtome sectioning. For Laboratory professional use only.

INSTRUCTIONS FOR USE

1. Chose a Cassette model specific to each tissue or biopsy size.
 - a) Smaller specimens can be contained with the use of biopsy foam pads, biopsy bags or wrapping if necessary.
 - b) Larger specimens may be trimmed to ensure a good fit for optimal processing.
2. Identify each cassette with the appropriate specimen/patient information during sampling; use identification methods that are resistant to histological solvents.
 - Use Laboratory marker pens for manual identification.
 - For printers use the Simport cassette selection guide in annex or our selection tool at <https://www.simport.com/en/cassette/1-cassette-selection-tool>
3. Place specimen in the cassette.
4. Close the lid. (a clicking sound will be heard and felt when a proper closure is obtained).
5. Process samples according to specific laboratory techniques.
6. Leave enough space between cassettes to allow adequate fluid transfer.
7. During the embedding of the tissue or biopsy sample, fill the back of the cassette thoroughly with paraffin to ensure proper support of the specimen in the microtome.
8. Avoid air pockets in hot paraffin when embedding before it solidifies.

Warning: Always check proper lid closure before placing cassette in processing equipment. Do not use if the hook is missing or damaged this presents a risk of opening and loss of biopsy. Do not overfill the cassette to prevent lid opening during processing. Safely dispose of infectious substances following applicable rules and regulations on waste disposals.

Please report any serious incident to Simport and to the competent authority of the Member State where you are established.

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