Simport® Scientific Cryovial® - Instructions for use

EN

All vials are made of polypropylene and sterilized by gamma radiation. They can withstand temperatures from +121°C down to -196°C, but should be stored only in the gas phase of liquid nitrogen. Some models contain silicone seals. T310-10A cap is made of polyethylene.

INSTRUCTIONS FOR USE

1. Fill the Cryovial® with a maximum of 90% of the tube's volume to permit sample expansion.
2. Verify that the Cryovial® and screw cap are completely dry before closing; liquid drops will impair the seal in the gas phase of liquid nitrogen.
3. Tighten the Cryovial® screw cap normally; over-tightening may squeeze the silicone seal out between the tube and cap and could increase the risk of leakage.
4. Avoid racks or stands that create partial insulating effect, the entire Cryovial® must freeze uniformly.
5. Move the Cryovial® periodically during freezing to reduce the risk of super cooling and subsequent quick ice formation.
6. Store the vials in the gas phase of liquid nitrogen. Immersion in the liquid may cause infiltration of LN2 into the vials which could result in leaking or shattering when brought to room temperature.
7. Use safety equipment: gloves, hoods, face shields, etc... when removing vials from liquid nitrogen containers.
8. To thaw most cells, warming should occur as rapidly as possible; place the Cryovial® in a water bath (37°C to 40°C), keeping the vials in constant motion.

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Simport® Scientific Inc., 2588, Bernard-Pilon, Beloeil (Québec) J3G 4S5 Canada