

S600 The MultiRack™

Made of acetal

A newly designed tube support that can be used all around the lab. The MultiRack™ is available in three models to accommodate a full range of laboratory test tubes and centrifuge tubes up to 30 mm in diameter. Although it is one of the most attractive racks available today, it offers all the advantages required by the modern laboratory. Made of highly resistant acetal, it will not shatter or stain in contact with most laboratory chemicals. No coating to worry about, which can chip, peel or rust in a water bath.

The MultiRack™ is compact, lightweight and stackable in order to save as much space as possible. This is why it is ideal for incubators, refrigerators, freezers, under lab hoods and on bench tops. Not only is it submersible but will also sink and maintain stability without tipping over.

The MultiRack™ is made of three-tiers to facilitate the insertion and stability of tubes. The base tier has rounded wells with drain holes. Convenient handles on each side of the rack will ensure a safe grip when carrying it around. Interlocking feet allow stacking. Series S600-13 will accommodate all tubes up to a diameter of 13 mm while series S600-16 will accept tubes up to 16 mm in diameter including 15 ml centrifuge tubes. Model S600-30 is perfect for accommodating up to 18 x 50 ml centrifuge tubes. Autoclavable at 121 °C for 20 minutes. Available in five attractive colors. Dimensions: 293 x 115 x 65 mm H (11 1/2 x 4 1/2 x 2 1/2 in. H)



S610 The MultiRack™ Jr.



Made of acetal

Also available is the MultiRack™ Jr. having the same features and benefits as the larger model but will hold half the number of tubes. A great acquisition when space is more limited.



6 models available to accommodate a full range of laboratory tubes up to 30 mm in diameter.

Rack Cat. #	Capacity	Rack Cat. #	Capacity	For tubes	Color	Qty/Pk	Qty/Cs
S600-13B	84	S610-13B	42	up to 13 mm	Blue	1	10
S600-13G	84	S610-13B	42	up to 13 mm	Green	1	10
S600-13L	84	S610-13B	42	up to 13 mm	Lilac	1	10
S600-13O	84	S610-13B	42	up to 13 mm	Orange	1	10
S600-13Y	84	S610-13B	42	up to 13 mm	Yellow	1	10

Rack Cat. #	Capacity	Rack Cat. #	Capacity	For tubes	Color	Qty/Pk	Qty/Cs
S600-16B	60	S610-16B	30	up to 16 mm	Blue	1	10
S600-16G	60	S610-16B	30	up to 16 mm	Green	1	10
S600-16L	60	S610-16B	30	up to 16 mm	Lilac	1	10
S600-16O	60	S610-16B	30	up to 16 mm	Orange	1	10
S600-16Y	60	S610-16B	30	up to 16 mm	Yellow	1	10

Rack Cat. #	Capacity	Rack Cat. #	Capacity	For tubes	Color	Qty/Pk	Qty/Cs
S600-30B	18	S610-30B	9	25 to 30 mm	Blue	1	10
S600-30G	18	S610-30B	9	25 to 30 mm	Green	1	10
S600-30L	18	S610-30B	9	25 to 30 mm	Lilac	1	10
S600-30O	18	S610-30B	9	25 to 30 mm	Orange	1	10
S600-30Y	18	S610-30B	9	25 to 30 mm	Yellow	1	10