

# Simport®

*Canadian and Family Owned Since 1975*



## **B721-1 SAMPLE TUBE 2.5 ML**

Made of polystyrene

This tube is used with the Roche Cobas Analyzer. It has perfect clarity for easy viewing of contents. This graduated false bottom tube has a capacity of 2.5 ml and will accept bar code labels.

Dimensions: 13 x 75 mm

Cat. #	Volume	Qty/Pk	Qty/Cs
B721-1	2.5 ml	1000	2000

[www.simport.com](http://www.simport.com)

# T404-3 Flange Plug Caps

for B721-1 Sample Tubes

Made of Polyethylene

These caps have two flexible flanges to ensure a leakproof seal. They will fit into test tubes and also in round cuvettes and centrifuges tubes.

Cat. #	For tubes with I.D.	Color	Qty/Pk
T404-3N	12 mm	Natural	1000
T404-3B	12 mm	Blue	1000
T404-3G	12 mm	Green	1000
T404-3R	12 mm	Red	1000
T404-3W	12 mm	White	1000
T404-3Y	12 mm	White	1000



## S600 The MultiRack™

Made of acetal polymer

This rack is perfectly suited for the B721-1 Sample Tube. Can also be used all around the lab for tubes up to 13 mm dia. Although it is one of the most attractive racks available today, it offers all the advantages required by the modern laboratory. Made of highly resistant acetal polymer, it will not shatter or stain in contact with most laboratory chemicals. No coating to worry about, which can chip, peel or rust in a water bath.

The MultiRack™ is compact, lightweight and stackable in order to save as much space as possible. This is why it is ideal for incubators, refrigerators, freezers, under lab hoods and on bench tops. Not only is it submersible but will also sink and maintain stability without tipping over.

The MultiRack™ is made of three-tiers to facilitate the insertion and stability of tubes. The base tier has rounded wells with drain holes. Convenient handles on each side of the rack will ensure a safe grip when carrying it around. Interlocking feet allow stacking. Available in five attractive colors. Dimensions: 293 x 115 x 65 mm H (11 1/2 x 4 1/2 x 2 1/2 in. H)

Two other formats are available for accepting larger tubes. Please call, or visit [www.simport.com](http://www.simport.com) for further details.



## S610 The MultiRack™ Jr.

Made of acetal

Also available is the MultiRack™ Jr. having the same features and benefits as the larger model but will hold half the number of tubes. A great acquisition when space is more limited.

Two other formats are available for accepting larger tubes. Please call, or visit [www.simport.com](http://www.simport.com) for further details.

Rack Cat. #	Capacity	Rack Cat. #	Capacity	For tubes	Color	Qty/Pk	Qty/Cs
S600-13B	84	S610-13B	42	up to 13 mm	Blue	1	10
S600-13G	84	S610-13G	42	up to 13 mm	Green	1	10
S600-13L	84	S610-13L	42	up to 13 mm	Lilac	1	10
S600-13O	84	S610-13O	42	up to 13 mm	Orange	1	10
S600-13Y	84	S610-13Y	42	up to 13 mm	Yellow	1	10

**Simport**<sup>®</sup>  
Scientific

2588 Bernard-Pilon  
Beloil, Qc J3G 4S5 Canada

Phone: (450) 464-1723

Fax: (450) 464-3394

Email: [info@simport.com](mailto:info@simport.com)

Website: [www.simport.com](http://www.simport.com)

Distributed by: