

Simport®

Canadian and Family Owned Since 1975

If you **TRULY** care about your sample,
let us help you **PROTECT** its integrity!



The
LOCKMAILER™

Microscope Slide Jar

www.simport.com

The LOCKMAILER™

Microscope Slide Jar

Jar made of Polypropylene
Closure made of High Density Polyethylene



At last a tamper evident multi purpose container for mailing, staining or storing microscope slides.

Constructed of extra-strong and clear polypropylene, it will hold up to 4 standard 3 x 1 in. or 75 x 25 mm slides vertically. Inside channels are slotted to keep slides safely separated. Perfect also for slide conveyors and specimen slide transport between the doctor's office and the lab.

It incorporates a unique tamper evident leakproof screw cap ensuring your peace of mind during transport or storage situations where someone might have manipulated your slides without your prior knowledge. Can also be used without the tamperproof locking mechanism. For color coding purposes, use a Capinsert™ (see T345 Series) that may be inserted on top of closure. Ten different colors are available.

The container is designed for maximum stability on a bench top while having an internal volume of only 12 ml.

Dimensions: 35 x 87 mm H (1 3/8 in. x 3 7/16 in. H)



M950-4MA

Can also be used as a slide staining jar

Cat. #	Color	Qty/Pk	Qty/Cs
M950-4MA	Magenta	100	500

T345

COLOR CODING CAPINSERT™

Made of polypropylene

The Capinsert™ is used to color code a LockMailer™ according to your specific needs. It is inserted on top of the closure and has a write-on frosted area for sample identification.



Cat. #	Color	Cat. #	Color	Qty/Pk	Qty/Cs
T345B	Blue	T345P	Pink	100	500
T345GY	Gray	T345R	Red	100	500
T345G	Green	T345V	Violet	100	500
T345L	Lilac	T345W	White	100	500
T345O	Orange	T345Y	Yellow	100	500
		T345AS	Assorted*	100	500

* Blue, lilac, red, yellow and white

HOW TO USE THE LOCKMAILER™



Push up the attached tab on side of vial.



Make sure it firmly clicks in place.



Screw on tamperproof cap all the way.



When opening the vial, the tamper evident ring will detach itself from the cap.



The LockMailer™ can also be used without the tamper evident feature.

Simport®
Scientific

2588 Bernard-Pilon
Beloeil, Qc J3G 4S5 Canada

Phone: (450) 464-1723

Fax: (450) 464-3394

Email: info@simport.com

Website: www.simport.com

Capinsert and LockMailer are trademarks of Simport® Scientific.

Distributed by: