



## S600 The MultiRack™

Made of acetal polymer

A newly designed tube support that can be used all around the lab. The MultiRack™ is available in three models to accommodate a full range of laboratory test tubes and centrifuge tubes up to 30 mm in diameter. Although it is one of the most attractive racks available today, it offers all the advantages required by the modern laboratory. Made of highly resistant acetal polymer, it will not shatter or stain in contact with most laboratory chemicals. No coating to worry about, which can chip, peel or rust in a water bath.

The MultiRack™ is compact, lightweight and stackable in order to save as much space as possible. This is why it is ideal for incubators, refrigerators, freezers, under lab hoods and on bench tops. Not only is it submersible but will also sink and maintain stability without tipping over.

The MultiRack™ is made of three-tiers to facilitate the insertion and stability of tubes. The base tier has rounded wells with drain holes. Convenient handles on each side of the rack will ensure a safe grip when carrying it around. Interlocking feet allow stacking. Series S600-13 will accommodate all tubes up to a diameter of 13 mm while series S600-16 will accept tubes up to 16 mm in diameter including 15 ml centrifuge tubes. Model S600-30 is perfect for accommodating up to 18 x 50 ml centrifuge tubes. Autoclavable at 121 °C for 20 minutes. Available in five attractive colors.


Dimensions: 293 x 115 x 65 mm H (11 1/2 x 4 1/2 x 2 1/2 in. H)

The MultiRack™ is available in 3 models to accommodate a full range of laboratory test tubes and centrifuge tubes up to 30 mm in diameter.

Cat. #	For tubes	Capacity	Color	Qty/Pk	Qty/Cs
S600-13B	up to 13 mm	84	Blue	1	10
S600-13G	up to 13 mm	84	Green	1	10
S600-13L	up to 13 mm	84	Lilac	1	10
S600-13O	up to 13 mm	84	Orange	1	10
S600-13Y	up to 13 mm	84	Yellow	1	10

Cat. #	For tubes	Capacity	Color	Qty/Pk	Qty/Cs
S600-16B	up to 16 mm	60	Blue	1	10
S600-16G	up to 16 mm	60	Green	1	10
S600-16L	up to 16 mm	60	Lilac	1	10
S600-16O	up to 16 mm	60	Orange	1	10
S600-16Y	up to 16 mm	60	Yellow	1	10

Cat. #	For tubes	Capacity	Color	Qty/Pk	Qty/Cs
S600-30B	25 to 30 mm	18	Blue	1	10
S600-30G	25 to 30 mm	18	Green	1	10
S600-30L	25 to 30 mm	18	Lilac	1	10
S600-30O	25 to 30 mm	18	Orange	1	10
S600-30Y	25 to 30 mm	18	Yellow	1	10



Also suitable for 15 ml centrifuge tubes and urine collection tubes.



This model is ideal for 50 ml centrifuge tubes.



## S510-500 The SecuRack™

Made of high impact polystyrene



Silicone tabs around each opening securely hold tubes in place.



This special 50-place rack will hold your 12 x 75 mm and 13 x 100 mm tubes securely in place thanks to silicone tabs surrounding the base of each tube while in the rack. This makes it very convenient to empty tube content before discarding them. Also great for holding tubes in rack securely in place when in a water bath. Each position is alpha numerically identified. Units can be anchored laterally to one another, thanks to two screws supplied with each rack. Dimensions: 250 x 128 x 50 mm H (9 3/4 x 5 x 2 in. H).

Cat. #	Color	Qty/Pk	Qty/Cs
S510-500	Orange	1	2