



SnapTwist™ Micrewtube®

Tube made of polypropylene
Cap made of polyethylene

The SnapTwist Micrew™ provides all the advantages of a modern microcentrifuge tube with screw cap but the SNAPTWIST® closure is a true time saver. The SnapTwist Micrew™ has a multitude of applications around your lab. It is ideal for freezer storage, boiling application, centrifugation etc, and will fit most standard microcentrifuge rotors.

The tubes can be securely sealed by simply capping the closures on. Removal of caps requires an easy 1/4 turn (twist). The ease with which these caps can be manipulated eliminates the danger of spillage associated with other push-on/pull-off caps. The deep internal lip of the cap fits snugly against the interior wall of the tube preventing the contents from coming in contact with the threads, thus reducing the chances of sample contamination. The quality of the sealing system is such that it is not necessary to tighten the closure with pressure to achieve a leakproof seal.

This series of tubes is not available with graduations but caps are supplied with attachment loops in order to prevent contamination and mix-up. Caps cannot be autoclaved since they are made of high density polyethylene. Conical bottom tubes can be centrifuged up to 20,000g. Temperature range: -90 °C to +100 °C



T345 Capinserts™ are available for a multitude of Simport products. For more details and colors available please refer to page 115.

SNAPTWIST™

Tubes and caps in the Micrewtube® Family are certified RNase, DNase, DNA and Pyrogen free



Two types of caps available. The one with attached loop helps avoid mix-ups and possible contamination.



All microcentrifuge tubes in the SnapTwist™ Family have a super clear highly polished surface for better viewing of contents.



Tubes

Cat. #	Description	Volume	Qty/Pk
T342-4T	Self-Standing	1.5 ml	1000
T342-5T	Conical Bottom	1.5 ml	1000
T342-6T	Self-Standing	1.8 ml	1000
T342-7T	Conical Bottom	1.8 ml	1000



Caps

Cat. #	Description	Color	Qty/Pk
T343NLS	Without Loop	Natural	1000
T343BLS	Without Loop	Blue	1000
T343GLS	Without Loop	Green	1000
T343LLS	Without Loop	Lilac	1000
T343RLS	Without Loop	Red	1000
T343YLS	Without Loop	Yellow	1000
T343WLS	Without Loop	White	1000



The vial can be securely sealed by simply snapping the cap on.



Cat. #	Description	Color	Qty/Pk
T343NLSL	With Loop	Natural	1000
T343BLSL	With Loop	Blue	1000
T343GLSL	With Loop	Green	1000
T343LLSL	With Loop	Lilac	1000
T343RLSL	With Loop	Red	1000
T343YLSL	With Loop	Yellow	1000
T343WLSL	With Loop	White	1000



Removal of cap requires an easy 1/4 turn (twist).