

C571 & C572

30 and 50 ml Sample Tubes

Tube made of polypropylene / Cap made of polyethylene


Conical bottom tubes with self-standing base graduated from 10 to 30 ml and 10 to 50 ml respectively in 5 ml increments. Chemically resistant and shatterproof, they are supplied with a leakproof screw cap, particularly important when transporting hazardous material. Ideal for transport and storage of urine, sputum and most liquids or particulate samples. Molded-in graduations make them easy to read. Available sterile or non sterile. The caps of Series C572 accept T345 CapInserts (see page 111 for details and colors). Functional temperature range: -90 °C to +121 °C. 30 ml tube dimensions: 25.3 mm dia x 111 mm length. 50 ml tube dimensions: 30 mm dia x 115 mm length.



T345 CapInserts™ are available for a multitude of Simport products. For more details and colors available please refer to page 111.



C572



C571

50 ml Sample Tube

Cat. #	Type	Volume (ml)	Color of cap	Qty/Bag	Qty/Cs
C571-1	Sterile	50	Green	25	500
C571-2	Non sterile	50	Yellow	100	500

30 ml Sample Tube

Cat. #	Type	Volume (ml)	Color of cap	Qty/Bag	Qty/Cs
C572-1	Sterile	30	White	25	500
C572-2	Non sterile	30	White	100	500

CONTAINERS

C580

Specimen Containers with Snap Cap

Container made of Polypropylene / Snap Cap made of Low Density Polyethylene

These disposable containers made of polypropylene, are ideal for collection and storage of fluids, powders, solids, pathology specimens and hazardous samples. Lids are tight fitting.

Cat. #	(ml/oz)	Opening (mm)	Height (mm)	Qty/Cs
C580-3	300/10	111	45	100
C580-4	500/17	111	75	100
C580-5	1000/34	111	140	100



The 2500 and 5500 ml sizes are provided with a plastic handle to facilitate handling



C581

Tamperproof Specimen Containers

Container made of polypropylene / Closure made of polyethylene

These disposable containers are ideal for collection, transport and storage of fluids, powders, solids, pathology specimens and hazardous samples. Tight-fitting lids prevent leaks and odors. They incorporate a tamperproof lid which is opened by first removing tab on side of lid, then by lifting the cover. They resist temperatures up to 121 °C. The two larger sizes are supplied with a plastic handle for easy carrying.

Cat. #	(ml/oz)	Opening (mm)	Height (mm)	Qty/Cs
C581-300	300/10	111	45	100
C581-500	500/17	111	75	100
C581-1L	1000/34	111	140	100
C581-2L	2500/80	150	168	25
C581-5L	5500/170	215	190	10