

M960

HistoTainer™ Prefilled Specimen Containers

Container made of polypropylene / Closure made of polyethylene

Especially designed for collection, transport and storage of histology specimens, Simport offers shatter resistant polypropylene containers, eliminating most problems of leakage and evaporation. Containers are uniquely stackable, shatter resistant and are manufactured from virgin, translucent polypropylene. The magenta lids are ribbed for easy opening when hands are wet or gloved while the jars are stackable for easy, safe storage and translucent to allow specimens to be viewed without opening. They are manufactured from virgin polyethylene with a unique integrated leak-resistant seal. Containers are available in the more rigid straight side format.

The uniqueness of the HistoTainer™ is that it incorporates an innovative tamper evident screw cap ensuring your peace of mind during transport or storage situations where someone might have manipulated the specimen without your prior knowledge. Can also be used without using the tamper evident locking mechanism. For color coding purposes, use a Capinsert™ (see T345 Series on page 48) on top of the closure. Ten different colors are available.

Both containers and caps are manufactured without the use of plasticisers or mold release agents. All material used in manufacture are free from latex. All containers are 95 kPa compliant. Available in many sizes from 20 to 120 ml.

The Simport HistoTainer™ is half filled with 10% Neutral Buffered Formalin as a fixative. 10% Neutral Buffered Formalin penetrates quickly, but fixes slowly. The Simport Formalin is enhanced by a buffering capacity that optimizing histological results by light microscopy and immunohistochemistry.

The Simport HistoTainer™ is highly recommended in the following laboratories:

- Microbiology
- Urinalysis
- Pathology
- Drug Testing
- Histology



The HistoTainer™ is built with a unique tamper evident locking mechanism used when sample integrity is of high importance.



Cat. #	Description	Volume	Qty/Pk	Qty/Cs
M960-20FMA	Histotainer Prefilled	20 ml	24	96
M960-40FMA	Histotainer Prefilled	40 ml	24	96
M960-60FMA	Histotainer Prefilled	60 ml	24	96
M960-90FMA	Histotainer Prefilled	90 ml	24	96
M960-120FMA	Histotainer Prefilled	120 ml	24	96



Packaged in sturdy cardboard boxes with handles for easy carrying

SPECIMEN CONTAINERS

Simport
Belsoil QC, Canada
10% NEUTRAL BUFFERED FORMALIN v/v
pH 6.8-7.2 @ 25°C

CAUTION: Contains Formaldehyde. Toxic by inhalation and if swallowed. Irritating to the eyes, respiratory system and skin. May cause sensitization by inhalation or skin contact. Risk of serious damage to eyes. Cancer Hazard-Repeated or prolonged exposure increases the risk.

Warning label has space for patient identification

Patient
 Date
 Specimen
 Doctor

T345

Color Coding CAPINSERT™



Made of polypropylene

T345 Capinserts™ are available for a multitude of Simport products. For more details and colors available please refer to page 111.

How to use the HistoTainer™



Remove Screw Cap.



Place sample in container.



Push up the attached tab on side of vial. Make sure it firmly clicks in place.



Screw on tamper evident cap completely.



When opening the vial, the tamper evident ring will detach itself from the cap.



The HistoTainer™ can also be used without the tamper evident feature.