

# Technical Data Sheet

|                                                    |                                   |
|----------------------------------------------------|-----------------------------------|
| <b>Product Name: CytoSep™</b>                      | <b>Catalogue No.: M984 Series</b> |
| <b>This document replaces any previous version</b> |                                   |

## 1. Product Description:

### ➤ Available Variations and/or Components:

- **M984-10FW** – Non-Sterile and Disposable Cytology Single-Funnel with Filter & Cap for 0.1 to 0.5 mL Sample Volumes
- **M984-20FW** – Non-Sterile and Disposable Cytology Double-Funnel with Filter & Cap for 0.05 to 0.25 mL Sample Volumes
- **M984-30M** – Non-Sterile and Disposable Single-Funnel with Filter & Cap for 0.5 to 6.0 mL Sample Volumes
- **M984FW10** – Non-Sterile and Disposable White filter paper for single funnel
- **M984FW20** – Non-Sterile and Disposable White filter paper for double funnel

### ➤ Compatible with:

- EpreDia™ (formerly Shandon) Cytospin™ 4 Cytocentrifuge and equivalent systems

## 2. Packaging:

- **M984-10FW:** 200 units per case
- **M984-20FW:** 200 units per case
- **M984-30M:** 25 units per case
- **M967-10FW:** 48 units per case
- **M967-20FW:** 48 units per case

## 3. Product Specifications:

### ➤ Material:

- **M984-10FW:** Polypropylene (PP)
- **M984-20FW:** Polypropylene (PP)
- **M984-30FW:** Polypropylene (PP)
- **M984FW10:** White, smooth surface, Cellulose
- **M984FW20:** White, smooth surface, Cellulose

## 4. Standards and Conformity:

- **FDA:** Resin conforms to FDA 21 CFR 177.1520  
Filter Papers conform to FDA 21 CFR 176.170 and 21 CFR 176.180
- **CE:** These Products are CE Marked
- **CH:** These Products are CH Marked
- **CONEG / RoHS:** Supplier certifies resin is in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Supplier certifies the resin is in conformity to REACH standards
- **EU Reg. 2019/1021 (POPs):** Supplier certifies the resin is in conformity with EU 2019/1021 (POP)
- **Toxic Substance Control Act:** Supplier certifies the resin is in conformity with US TSCA
- **PFAS:** Supplier certifies the resin is PFAS-Free
- **BSE/TSE:** Supplier certifies the resin is BSE/TSE-Free
- **Bisphenols:** Supplier certifies the resin is free from Bisphenols & their compounds
- **Phthalates:** Supplier certifies the resin is Phthalate Free
- **LATEX:** Supplier certifies the resin is Latex Free
- **Conflict Minerals Reporting:** Supplier certifies the resin is in conformity with CMR standards

# Technical Data Sheet

## 5. Quality Assurance:

- No presence of contamination in plastic
- Assembly and fit
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any damage to the packaging

## 8. Recommended Use:

- Follow chemical resistance chart recommendations
- Check cap closure when using biohazard material and / or chemical reagents

For use in EpreDia™ Cytospin™ 4 Cyto centrifuge automated equipment, please follow the equipment manufacturer's instructions for use

## **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*

## **Disclaimer:**

*The information contained in this document shall not constitute a representation or warranty of the suitability for any specific purpose. To the best of our knowledge, the information contained is reliable and true based on our internal and external testing as well as statements supplied to us by third party suppliers.*