

Product Name: SecurTainer™ I

Catalogue No.: C575-60MA

This document replaces any previous version

1. Product Description:

- **60ml Specimen Container:** Disposable, non-sterile with tamper evident leakproof screw cap

2. Packaging:

- **Case:** 5 packs of 100 units / 500 units per case

3. Product Specifications:

- **Material:**
 - Container: Polypropylene
 - Cap: Polyethylene
 - Cap Color: magenta
- Caps packaged separately in bags of 100
- Gamma radiation resistant
- Temperature range: -90°C to +100°C
- Cap configuration allows insertion of a Capinsert™ (T345 Series) for color-coding

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resins conform to FDA 21 CFR 177.1520
- **RoHS** Plastics and colorants are in conformity with RoHS standards for heavy metals
- **REACH (SVHC):** Plastics are in conformity to REACH standards
- **LATEX:** Materials are Latex Free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension and volume measurements
- Closure verification
- Leak proof testing in vacuum at 71.3cm Hg

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*