

Product Name: Petri Dish

Catalogue No.: D210-17

This document replaces any previous version

1. Product Description:

- **20ml Contact Plate:** Disposable, sterile petri dish with alphanumeric grid on lid

2. Packaging:

- **Case:** 25 sleeves of 20 units / 500 units per case

3. Product Specifications:

- **Material:**
 - Dish: Polystyrene
- **Actual dimension:** 60 mm x 15 mm
- **Sterility:** All lots are manufactured using aseptic process which has been validated in accordance with USP <71>
- **Temperature range:** -80°C to +80°C
- **Non autoclavable**
- **10 mm² alphanumeric grid**

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640
- **CONEG:** Plastic and colorant are in conformity with CONEG standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free
-

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging which could compromise the product sterility

8. Recommended Use:

- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*