

Edition Date: Revision Date : Revised by:

Product Name: Antistatic Weighing Dish Catalogue No.: D251-1 This document replaces any previous version

1. Product Description:

12 ml Pour Boat Dish: Disposable, non-sterile, antistatic, flat bottom with a smooth surface to facilitate dispensing of powdered and liquid samples

2. Packaging:

Case: 8 packages of 250 units / 2 000 units per case

3. Product Specifications:

- > Material:
- Dish: Antistatic Polystyrene
- Pour Boat Configuration
- Temperature range -80°C to +80°C
- Not autoclavable

4. Standards and Conformity:

➢ ISO 2859-1:	Sampling and inspection procedures
---------------	------------------------------------

- FDA: Resin conforms to FDA 21 CFR 177.1640
- **REACH (SVHC):** Plastic is in conformity to REACH standards

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- > The lot number can be found in one or all of these locations:
 - 1. On exterior case label
 - 2. On label inserted inside the master case
 - 3. On the inner bag



Edition Date: Revision Date : Revised by: 2007-09-14 2020-08-28 Annette Roy, Regulatory

7. Storage Conditions:

- > Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- > Follow chemical resistance chart recommendations
- > For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

Materials and information contained in this document are **Simport Scientific**'s copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.