

Product Name: **Antistatic Weighing Dish**

Catalogue No.: **D251-3**

This document replaces any previous version

1. Product Description:

- **270 ml Pour Boat Dish:** Disposable, non-sterile, antistatic, flat bottom, smooth surface and white pour boat to facilitate dispensing of powdered and liquid samples

2. Packaging:

- **Case:** 250 units per case

3. Product Specifications:

- **Material:**
 - Dish: Antistatic Polystyrene
- **Temperature range** -80°C to +80°C
- **Not autoclavable**

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640
- **REACH (SVHC):** Plastic is in conformity to REACH standards

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag



Simport
Since 1975 *Scientific inc.*

Technical Data Sheet

Edition Date: 2007-09-14
Revision Date : 2020-08-28
Revised by: Annette Roy, Regulatory

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*