

**Product Name: UNISSETTE**

**Catalogue No.: M405 Series**

This document replaces any previous version

## 1. Product Description:

- **Tissue Cassette with Attached Lid:** Non-Sterile, disposable cassette compatible with Primera Printer
- **Available Colors:**
  - White, Pink, Green, Yellow, Blue, Peach, Tan, Gray, Lilac, Orange and Aqua
  - Fluorescent Pink, Green and Yellow also available, minimum order quantity required

## 2. Packaging:

- **Case:** 3 inner packs of 500 units / 1 500 units per case

## 3. Product Specifications:

- **Material:**
  - Acetal
- **Writing surface:** Smooth surface with 35° angle
- **Compatible with Primera cassette printers**

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.2470
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free
- **CE:** Product is CE marked

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Lid closure verification

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

## 8. Recommended Use:

- Verify proper lid closure of the cassette before using
- Avoid picric acid, prolonged use in this acid will cause degradation
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*