

**Product Name:** **Metal lid**

**Catalogue No.:** **M481**

This document replaces any previous version

## 1. Product Description:

- **Metal lid:** Non-sterile lid for M480 cassettes

## 2. Packaging:

- **Case:** 25 units per package

## 3. Product Specifications:

- **Material:**  
Stainless Steel

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

## 5. Quality Assurance:

- No presence of contamination on stainless steel
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions

- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

## 8. Recommended Use:

- Accessory to be used in combination with cassettes for embedding procedures

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*