

**Product Name:** **Modular Storage Drawer**

**Catalogue No.:** **M495-6**

This document replaces any previous version

## 1. Product Description:

- **Modular Storage Drawer:** Non-sterile, capacity of 165 embedding rings or 250 cassettes per drawer

## 2. Packaging:

- **Case:** 6 units per case

## 3. Product Specifications:

- **Material:**
  - Drawer: High Impact Polystyrene
- 6 identification labels included per package
- Stackable with interlocking ridges on top and bottom
- Gamma radiation resistant
- Temperature range: -80°C to +80°C

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals

## 5. Quality Assurance:

- No presence of contamination in plastic
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

**7. Storage Conditions:**

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

**8. Recommended Use:**

- Follow chemical resistance chart recommendations

***OWNERSHIP OF MATERIALS:***

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*