

**Product Name:** SlideFolder™

**Catalogue No.:** M750-20 Series

This document replaces any previous version

## 1. Product Description:

- **Slide Storage System:** Non-sterile, stackable slide holder
- **Available Colors:**
  - Blue, Green, Gray, White, Yellow and assorted (2 of each color)

## 2. Packaging:

- **Case:** 10 units per case

## 3. Product Specifications:

- **Material:**
  - Base: High-Impact Polystyrene
  - Cover: Polystyrene
- Can hold up to 20 microscope slides of 75mm x 25mm (3 x 1 in.)

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

## 8. Recommended Use:

- None

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*