

**Product Name:** Cryosette™

**Catalogue No.:** M9560

This document replaces any previous version

## 1. Product Description:

- **Tissue Storage Container with Lid:** Disposable and non-sterile
- **Available Colors:**
  - Orange container
  - White lid

## 2. Packaging:

- **Case:** 250 containers with 5 bags of 50 lids per case / 250 units per case

## 3. Product Specifications:

- **Material:**
  - Acetal container
  - High Density Polyethylene lid
- Slanted Writing surface
- Temperature range: -80°C to +30°C

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.2470 and CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- closure verification

**6. Traceability:**

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

**7. Storage Conditions:**

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging

**8. Recommended Use:**

- Verify proper lid closure when using biohazard material
- Avoid picric acid, prolonged use in this acid will cause degradation
- Follow chemical resistance chart recommendations

**OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*