Edition Date:

Revision Date:

Revised by: Annette Roy, Regulatory

2021-12-17

Product Name: CytoSep™ Catalogue No.: M963-1

This document replaces any previous version

- 1. **Product Description:**
 - > Fluid Chamber: Disposable cytology Funnel, 1 ml fluid chamber
- 2. Packaging:
 - Case: 4 packages of 50 units / 200 units per case
- 3. **Product Specifications:**

Compatible with the Sakura Cyto-Tek® Cytocentrifuge:

- Material:
 - Polypropylene
- 4. **Standards and Conformity:**

> ISO 2859-1: Sampling and inspection procedures

Material is Latex Free > LATEX: **≻** CE: Product is CE marked

- 5. **Quality Assurance:**
 - Clear, no presence of contamination in plastic
 - Visual attributes
- 6. Traceability:
 - **➤** Lot No. Composition: 8 or 9 digits
 - > The lot number can be found in one or all of these locations:
 - 1. On exterior case label
 - 2. On label inserted inside the master case
 - 3. On the inner bag

Page 1 / 2 Issue Date: 2017-06-15

Edition Date: Revision Date: 2021-12-17

Revised by:

Annette Roy, Regulatory

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any possible contamination

8. Recommended Use:

For use in **Sakura Cyto-Tek® Cytocentrifuge** automated equipment, follow the equipment manufacturer's instructions for use

OWNERSHIP OF MATERIALS:

Materials and information contained in this document are **Simport Scientific**'s copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.

Document Template: FRM-106 Rev01 Technical Data Sheet
Page 2 / 2

Issue Date: 2017-06-15