

Product Name: **CytoSep™**

Catalogue No.: **M964-20FW1**

This document replaces any previous version

1. Product Description:

- **Double Cytology Funnel with Cap:** Disposable, pre-assembled with white paper filter

2. Packaging:

- **Case:** 100 individually packaged per inner pack / 500 units per case

3. Product Specifications:

Compatible with the *Shandon Cytopspin®* Cytocentrifuge:

- **Material:**
 - Funnel: Polypropylene
 - Cap: Polypropylene
 - Filter: White Paper

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **LATEX:** Material is Latex Free
- **CE:** Product is CE marked

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination

8. Recommended Use:

- For use in **Shandon Cytospin® Cyto centrifuge** automated equipment, follow the equipment manufacturer's instructions for use
- The double cytology funnel allows for two samples to be run simultaneously on a single slide

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*