

Product Name: OneHand™ Rack

Catalogue No.: T360-50

This document replaces any previous version

1. Product Description:

- **Rack:** Non-sterile, 50 position, stackable rack with serrated holes for tubes with ribs and locking base for one-hand operations with skirted bottom Simport Micrewtube® 0.5 and 2 ml

- **Available colors:**
 - Blue, Green, Lilac, Orange and Yellow

- **Compatible with:**
 - 0.5 to 2 ml Micrewtubes®

2. Packaging:

- **Case:** 10 units per case

3. Product Specifications:

- **Material:**
 - Acetal
- Temperature range: -80°C to +30°C
- Not Autoclavable
- Total capacity: 50 Microtubes

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions;

8. Recommended Use:

Rack for holding tubes

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*