

Edition Date: Revision Date : Revised by:

Product Name: Septum Screw Caps Catalogue No.: T504AQX

This document replaces any previous version

- **1.** Product Description:
 - Septum Screw Caps: Disposable, non-sterile cap with pierceable septum for external threaded sample tubes
 - > Compatible with:
 - T501 Sample Tubes

2. Packaging:

Case: 250 units per package

3. Product Specifications:

- > Material:
 - Cap: Polypropylene
 - Septum top: Teflon
 - Septum inner liner: Silicone
- Temperature range: -80°C to +121°C
- > Autoclavable at 121°C, for up to 30 minutes

4. Standards and Conformity:

- ISO 2859-1: Sampling and inspection procedures
- ➢ LATEX: Material is Latex Free
- **BSE / TSE:** Material is BSE / TSE Free

Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements
- Closure verification



Edition Date: Revision Date : Revised by: 2013-08-21 2021-10-26 Annette Roy, Regulatory

5. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- > The lot number can be found in one or all of these locations:
 - 1. On exterior case label
 - 2. On label inserted inside the master case
 - 3. On the inner bag

6. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging

7. Recommended Use:

- > Follow chemical resistance chart recommendations
- > For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

Materials and information contained in this document are **Simport Scientific**'s copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.