

Product Name: Uritainer®

Catalogue No.: B360-35

This document replaces any previous version

1. Product Description:

- **3.5 Liters, 24-Hr Urine Container:** Disposable, Non-sterile assembled with screw cap with washer liner.

2. Packaging:

- **Case:** 24 units per case

3. Product Specifications:

- **Material:**
 - Bottle: Polyethylene
 - Cap: Polypropylene
 - Valve: Polyethylene
- Gamma radiation resistant
- Temperature range: -80°C to +100°C.

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG /RoHS standards for heavy metals
- **LATEX:** Material is Latex free

5. Quality Assurance:

- No presence of contamination in plastic;
- Visual attributes;
- Closure and valve verification.

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case;
 3. On the inner bag.

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*