

Product Name: Tricorn™ Beaker

Catalogue No.: B700 series

This document replaces any previous version

1. Product Description:

- **Beakers:** Disposable, non-sterile

Product code	Volume (ml)	Graduated increments (ml)
B700-50	20	5
B700-100	100	10
B700-250	250	10
B700-400	400	20
B700-800	800	50
B700-1L	1000	50

2. Packaging:

- **Case:** 4 bags of 25 units / 100 units per case

3. Product Specifications:

- **Material:** Polypropylene
- **Temperature range:** -80°C to +121°C
- **Autoclave at +121°C,** for up to 30 minutes

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** plastics conform to FDA 21 CFR 177-1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals

5. Quality Assurance:

- No presence of contamination in plastic
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions

8. Recommended Use:

- Follow chemical resistance chart recommendations of Polypropylene

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*