

**Product Name: Tricorn™ Beaker**

**Catalogue No.: B700-1L**

**This document replaces any previous version**

## 1. Product Description:

- **1000 ml Beaker:** Disposable, non-sterile

## 2. Packaging:

- **Case:** 4 bags of 25 units / 100 units per case

## 3. Product Specifications:

- **Material:**
  - Beaker: Polypropylene
- **Volume:** 1000 ml
- **Graduated:** 50 ml increments
- **Temperature range:** -80°C to +121°C
- **Autoclave at +121°C,** for up to 30 minutes

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals

## 5. Quality Assurance:

- No presence of contamination in plastic
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions

## 8. Recommended Use:

- Follow chemical resistance chart recommendations of Polypropylene

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*