

Edition Date:	2022-04-22
Edited by:	Annette Roy, Regulatory and Technical Support
Revision Date :	2022-07-14
Revised by:	Christine Laforge RT, Technical Representative

Product Name: **MicrewTube™ Cap**

Catalogue No.: **T366LS Series**

This document replaces any previous version

1. Product Description:

- **MicrewTube™ Cap:** Non-sterile, cap with lip seal for MicrewTube™ Tubes sold separately.
- **Available colors:**
 - Natural, Amber, Blue, Green, Red, Yellow and White
- **Compatible with:**
 - T365 MacrewTube™ Series

2. Packaging:

- **Case:** 2 packs of 100 / 200 units per package

3. Product Specifications:

- **Material:** Polyethylene cap
- **Gamma sterilization resistant**
- **Temperature range:** -196°C to +100°C

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **CONEG:** Plastics and colorants are in conformity with CONEG standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex free
- **BSE / TSE:** Material is BSE / TSE free



Technical Data Sheet

Edition Date:	2022-04-22
Edited by:	Annette Roy, Regulatory and Technical Support
Revision Date :	2022-07-14
Revised by:	Christine Laforge RT, Technical Representative

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements
- Closure verification

6. Traceability:

- **Lot No. Composition:** 8 to 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*