

## TECHNICAL DATA SHEET

### SCL21520

#### **PRODUCT CODE**

SCL21520

#### **PACKAGING**

Case of 500 (2 boxes of 250 bags)

#### **SHELF LIFE**

Five years sterility guaranteed

#### **PRODUCT DESCRIPTION**

*SECURE-T sterile sampling bags provide a sterile, secure and flexible receptacle that ensures reliable analysis results.*

*The bags can accommodate various types of solid, semi-solid and liquid samples. Easy to use and with multiple possible uses.*

#### **PRODUCT SPECIFICATIONS**

Closure type	<i>Tear-Off protection strip</i>
Closure material	<i>N/A</i>
Printed	<i>Clear – Non printed material</i>
MSDS	<i>Not applicable for this product</i>
Sterilisation	<i>Sterility is assured due to extrusion temperatures exceeding 220°C which eliminates any present bacteria. During manufacturing all bags are sealed at both top and bottom in a continuous format before being cut into individual bags. This method assures that the bag's interior will never come in contact with its surroundings. Periodic validation tests, based on USP 40 – Sterility testing, are performed in our facility to guaranty the sterility of our processes.</i>
RNase/DNase/Pyrogen	<i>Products are free of RNase, DNase and pyrogens</i>
Volume	<i>405 oz (12,000mL)</i>
Dimensions	<i>15 in x 20 in (381 mm x 508 mm)</i>
Thickness of material	<i>2 mil (51 micron)</i>
FDA regulation	<i>All TWIRL'EM sterile sampling bags comply to the applicable regulations of the Food and Drug Administration (21 CFR 177.1520 (c)3.1a, (c)3.2a and 178.2010.</i>
EU regulation	<i>All TWIRL'EM sterile sampling bags comply where applicable to EU No 10/2011 regulation for use in contact with food, subject to certain types of migration limitations and labeling requirements.</i>
CFIA	<i>All TWIRL'EM sterile sampling bags are also a solution which could be used in the CFIA Prevention Control Plan inspired by HACCP principles, a Canadian federal initiative from the Safe Food for Canadians Regulations (SFCR).</i>

T : 450-649-7343

F : 450-649-3113

www.labplas.com

service@labplas.com

Update : 14/4/2020