

**Product Name:** **Sample Cup**

**Catalogue No.:** **B720-16**

This document replaces any previous version

## 1. Product Description:

- **1 ml Sample Cup:** Disposable, Non-sterile, 16.7 mm diameter sample cup, for use with SeraNest™.

## 2. Packaging:

- **Case:** 5 bag of 1 000 units / 5000 per case

## 3. Product Specifications:

- **Material:** Polystyrene
- **Color:** Natural
- **Temperature range for tube only:** -80°C to +80°C
- **Not Autoclavable**

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex free

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination

## 8. Recommended Use:

- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*