

Product Name: **SlideTray™**

Catalogue No.: **M755-20W**

1. Product Description:

- **The SlideTray™ is microscope slide holder.**
 - Easy to insert and remove slides
 - Holds up to 20 microscope slides
 - Stackable

 - **Available Colors:**
 - White

 - **Compatible with:**
 - Fits standard 3 x 1 inch microscope slides

2. Packaging:

- **Case:** 10 / units per case

3. Product Specifications:

- **Material:**
 - High Impact Polystyrene
- Temperature range: -80°C to +80°C
- Not autoclavable

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Closure verification

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging or product

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*