

Product Name: Flange Plug Caps

Catalogue No.: T404-10 Series

This document replaces any previous version

1. Product Description:

- **Flange Plug Cap for 16 mm Culture Tubes:** Non sterile, disposable cap
- **Available Colors:**
 - Natural, Blue, Green, Red, White and Yellow

2. Packaging:

- **Case:** 1 000 units per package

3. Product Specifications:

- Material: Polyethylene
- Temperature range: -80°C to +80°C
- Cap configuration has two flexible flanges
- Fits with 16 mm test tubes, round cuvettes, and centrifuge tubes

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **CONEG:** Plastics and colorants are in conformity with CONEG standards for heavy metals
- **Latex:** Material is Latex free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements
- Closure verification

6. Traceability:

- **Lot No. Composition:** 8 to 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*