

Product Name: Dropette®

Catalogue No.: P200-301S

This document replaces any previous version

1. Product Description:

- **Transfer Pipet:** Sterile, disposable pipet

2. Packaging:

- **Case:** packaged individually / 400 per inner pack / 4000 units per case

3. Product Specifications:

- **Material:**
 - Low Density Polyethylene
- Graduated increments from 0.25ml to 1ml
- Bulb draw: 2.1ml
- Total volume: 3ml
- Length: 13.8cm
- Temperature range : -80°C to +80°C
- Gamma radiation sterilized at a SAL of 10^{-3} ; specified dose between 6.5 kGy and 13.5 kGy

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **Latex:** Material is Latex free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension and volume measurements

6. Traceability:

- **No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any damage to the packaging

8. Recommended Use:

- Be cautious when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*