

Scientific inc. Technical Data Sheet

Edition Date: 2023-02-02 **Revision Date:** None

Revised by: Annette Roy, Regulatory and

Technical Support

Product Name: SuperFLO Catalogue No.: M430-T

Routine, Taped Cassette with attached lid

Material: Acetal

Case: 25 stacks of 40 taped cassettes / 1000 total units per case

Description: Disposable, Non-sterile

With attached lid, a 45° angle and compatible with cassette printers

Temperature range: -80°C to +30°C ISO 2859-1: Sampling and inspection

FDA: Plastics conform to FDA 21 CFR 177.2470 **CONEG / RoHS:** Plastics and colorants comply

LATEX: Material is Latex free CE: Product is CE marked

Recommended Use: For histological sample containment and processing

Storage Conditions: Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging.

OWNERSHIP OF MATERIALS: Materials and information contained in this document are **Simport Scientific**'s copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.

Document Template: FRM-106 Rev02 Technical Data Sheet Page 1 / 1

Issue Date: 2021-08-19