

Technical Data Sheet

Product Name: SuperFLO

Catalogue No.: M531

Biopsy Cassette with attached Lid

Material: Acetal

Case: 500 units per pack, 1500 per case

Description: Disposable, Non-sterile
With attached lid, a 45° angle and a surface ideal for handwriting

Temperature range: -80°C to +30°C

ISO 2859-1: Sampling and inspection

FDA: Plastics conform to FDA 21 CFR 177.2470

CONEG / RoHS: Plastics and colorants comply

LATEX: Material is Latex free

CE: Product is CE marked

Recommended Use: For histological sample containment and processing

Storage Conditions: Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging.

OWNERSHIP OF MATERIALS: *Materials and information contained in this document are Simport Scientific's copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*