

Scientific inc. Technical Data Sheet

Edition Date: 2023-02-02 Revision Date: None

Revised by: Annette Roy, Regulatory and

Technical Support

Product Name: SuperFLO Catalogue No.: M533-SL

Routine Cassette with separate hinged Lid

Material: Acetal

Case: 10 sleeves of 75 cassettes with 10 bags of 75 lids / 750 Total units

Description: Disposable, Non-sterile

With separate lid, a 45° angle and compatible with Thermo Fisher cassette printers

Temperature range: -80°C to +30°C **ISO 2859-1**: Sampling and inspection

FDA: Plastics conform to FDA 21 CFR 177.2470 **CONEG / RoHS:** Plastics and colorants comply

LATEX: Material is Latex free **CE:** Product is CE marked

Recommended Use: For histological sample containment and processing

Storage Conditions: Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging.

OWNERSHIP OF MATERIALS: Materials and information contained in this document are **Simport Scientific**'s copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.

Document Template: FRM-106 Rev02 Technical Data Sheet
Page 1 / 1

Issue Date: 2021-08-19