

Technical Data Sheet

Product Name: SuperFLO

Catalogue No.: M534-SL

Biopsy Cassette with separate hinged Lid

Material: Acetal

Case: 10 sleeves of 75 cassettes with 10 bags of 75 lids / 750 Total units

Description: Disposable, Non-sterile
With separate lid, a 45° angle and compatible with Thermo Fisher cassette printers

Temperature range: -80°C to +30°C

ISO 2859-1: Sampling and inspection

FDA: Plastics conform to FDA 21 CFR 177.2470

CONEG / RoHS: Plastics and colorants comply

LATEX: Material is Latex free

CE: Product is CE marked

Recommended Use: For histological sample containment and processing

Storage Conditions: Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging.

OWNERSHIP OF MATERIALS: *Materials and information contained in this document are Simport Scientific's copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*