

**1. Product Name: ANTISTATIC  
HEXAGONAL WEIGHING DISHES**

**Catalogue No.: D252-2**

## Weighing Dishes

**Material:** Antistatic polystyrene

**Case:** 250 units per pack, 500 per case

**Description:** Weighing Dishes. Volume: 58 ml

**Temperature range:** Up to 80°C

**ISO 2859-1:** Sampling and inspection

**FDA:** Plastics comply to FDA 21 CFR 177.1640

**REACH:** Plastic comply

**RoHS:** Plastic comply

**CONEG:** Comply

**LATEX:** Material is Latex free

**BSE/TSE FREE**

**Recommended Use:** They are ideal for many applications, such as weighing and dispensing.

**Storage Conditions:** Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging.

**OWNERSHIP OF MATERIALS:** *Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*