

Product Name: Specimen Container

Catalogue No.: C570-12

This document replaces any previous version

1. Product Description:

- **4 oz Specimen container:** Disposable, Non-sterile Specimen container with separate lids.

2. Packaging:

- **Case:** 10 packs of 100 units / 500 units per case

3. Product Specifications:

- **Material:**
 - Tube: Polypropylene
 - Cap: Polyethylene
- **Temperature range:** -90°C to +100°C
- **Tubes have molded-in graduations up to 128 ml**

4. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Closure verification

5. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

6. Storage Conditions:

- Store at room temperature in normal warehouse conditions

7. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*