

Technical Data Sheet

Product Name: **Histosette®**

Catalogue No.: **M393T Series**

This document replaces any previous version

1. Product Description:

➤ Biopsy Cassette with Attached Lid:

- Non-sterile, disposable cassette

➤ Available Colors:

- Standard Cassette (M393-XT) Colour options include:
 - M393-2T – White
 - M393-3T – Pink
 - M393-4T – Green
 - M393-5T – Yellow
 - M393-6T – Blue
 - M393-7T – Peach
 - M393-8T – Tan
 - M393-9T – Gray
 - M393-10T – Lilac
 - M393-11T – Orange
 - M393-12T – Aqua
- Fluorescent Pink, Green and Yellow also available, minimum order quantity required

2. Packaging:

- **Case:** **M393-XT:** 50 Quickload™ Stacks of 40 cassettes / 2 000 units per case

3. Product Specifications:

➤ Material:

- Acetal

➤ Writing Surface angle: 45° with specially designed surface for cassette printers

4. Standards and Conformity:

- **FDA:** Resin conforms to FDA 21 CFR 177.2470
- **CE:** These Products are CE Marked
- **CH:** These Products are CH Marked
- **CONEG / RoHS:** Supplier certifies resin is in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Supplier certifies the resin is in conformity to REACH standards
- **EU Reg. 2019/1021 (POPs):** Supplier certifies the resin is in conformity with EU 2019/1021 (POP)
- **Toxic Substance Control Act:** Supplier certifies the resin is in conformity with US TSCA
- **PFAS:** Supplier certifies the resin is PFAS-Free
- **BSE/TSE:** Supplier certifies the resin is BSE/TSE-Free
- **Bisphenols:** Supplier certifies the resin is free from Bisphenols & their compounds
- **Phthalates:** Supplier certifies the resin is Phthalate Free
- **LATEX:** Supplier certifies the resin is Latex Free
- **Conflict Minerals Reporting:** Supplier certifies the resin is in conformity with CMR standards

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Closure verification

Technical Data Sheet

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*

Disclaimer:

The information contained in this document shall not constitute a representation or warranty of the suitability for any specific purpose. To the best of our knowledge, the information contained is reliable and true based on our internal and external testing as well as statements supplied to us by third party suppliers.