

**Product Name: Disposable Funnel**

**Catalogue No.: F490-2**

This document replaces any previous version

## 1. Product Description:

- **Funnel:** Disposable and non-sterile

## 2. Packaging:

- **Case:** 4 stacks of 25 units / 100 units per case

## 3. Product Specifications:

- **Material:** Polystyrene
- **Top Inner Diameter:** 64 mm
- **Total Length:** 75.7 mm
- **Use with Paper Diameter:** 12.5 cm
- **Tapered at a 60° angle with inside fluting**

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** plastics conform to FDA 21 CFR 177-1640
- **CONEG:** Plastics and colorants are in conformity with CONEG standards for heavy metals

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any possible contamination
- Protect from any damage to the packaging

## 8. Recommended Use:

- Careful for spills when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*