**Edition Date:** 2017-08-28 Revision Date: N/A

Revised by: Annette Roy, QA Coordinator

Product Name: Base Mold Catalogue No.: M474

This document replaces any previous version

- **Product Description:** 1.
  - > Base Mold: Non-sterile, to fit with most cassettes and embedding rings
- 2. Packaging:
  - Case: 12 units per package
- 3. **Product Specifications:** 
  - Material:

Stainless Steel

Offered in 5 different sizes to accommodate each Biopsy or tissue size

M474-1 7 x 7 x 5 M474-2 15 x 15 x 5 M474-3 24 x 24 x 5 M474-4 30 x 24 x 5 M474-5 37 x 24 x 5

- 4. **Standards and Conformity:** 
  - > ISO 2859-1: Sampling and inspection procedures
- 5. **Quality Assurance:** 
  - ➤ No presence of contamination on stainless steel
  - Visual attributes
  - Dimension measurements
- 6. Traceability:
  - **Lot No. Composition:** 8 or 9 digits
  - > The lot number can be found in one or all of these locations:
    - 1. On exterior case label
    - 2. On label inserted inside the master case
    - 3. On the inner bag

Page 1 / 2 Issue Date: 2017-06-15

**Edition Date:** 2017-08-28 **Revision Date:** N/A

**Revised by:** Annette Roy, QA Coordinator

## 7. Storage Conditions:

- > Store at room temperature in normal warehouse conditions
- > Avoid temperature variations and humidity
- > Protect from any possible contamination
- Protect from any damage to the packaging

## 8. Recommended Use:

> Accessory to be used in combination with cassettes for embedding procedures

## **OWNERSHIP OF MATERIALS:**

Materials and information contained in this document are **Simport Scientific**'s copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.

Document Template: FRM-106 Rev01 Technical Data Sheet
Page 2 / 2

Issue Date: 2017-06-15