**Edition Date:** 2008-04-03 **Revision Date:** 2020-08-26

**Revised by:** Annette Roy, Regulatory

Product Name: Cryovial® Catalogue No.: T311-5

This document replaces any previous version

## 1. Product Description:

> 5 ml Cryogenic Vial with Cap: Sterile, disposable, , Round bottom, internal threaded design

assembled with washer seal cap

2. Packaging:

➤ Case: 10 bags of 100 units / 1 000 units per case

# 3. Product Specifications:

➤ Material:

Tube: PolypropyleneCap: Polypropylene

• Washer Seal: White Silicone

Certified RNase, DNase, Pyrogen and DNA Free

➤ Gamma radiation sterilized at a SAL of 10<sup>-3</sup>; specified dose between 6.5 kGy and 13.5 kGy

> Temperature range: -196°C to +121°C

➤ Autoclavable at +121°C for up to 30 minutes

> Tubes have printed graduations and marking area

➤ Cap configuration allows insertion of a Capinsert<sup>TM</sup> (T312 Series) for color-coding

## 4. Standards and Conformity:

➢ ISO 2859-1: Sampling and inspection procedures➢ FDA: Resin conforms to FDA 21 CFR 177.1520

> USP: Resin conforms to USP Class VI

> CONEG / RoHS: Plastics and colorants are in conformity with

CONEG / RoHS standards for heavy metals

> REACH (SVHC): Plastic is in conformity to REACH standards

➤ LATEX: Material is Latex Free
➤ BSE / TSE: Material is BSE / TSE Free
➤ CE: Product is CE marked

Document Template: FRM-106 Rev01 Technical Data Sheet

Issue Date: 2017-06-15

**Edition Date:** 2008-04-03 **Revision Date:** 2020-08-26

Revised by: Annette Roy, Regulatory

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Volume measurements
- Closure verification
- Leak proof testing in vacuum at 71.3cm Hg
- Gas phase of Liquid Nitrogen resistance

# 6. Traceability:

**Lot No. Composition:** 8 or 9 digits

## > The lot number can be found in one or all of these locations:

- 1. On exterior case label
- 2. On label inserted inside the master case
- 3. On the inner bag

# 7. Storage Conditions:

- > Store at room temperature in normal warehouse conditions
- > Avoid temperature variations and humidity
- > Protect from any possible contamination
- > Protect from any damage to the packaging which could compromise the product sterility

#### 8. Recommended Use:

- > Verify proper cap closure when using biohazard material and / or chemical reagents
- > Follow chemical resistance chart recommendations
- > For use in automated equipment, follow the equipment manufacturer's instructions
- > Should be used only in the gas phase of Liquid Nitrogen

#### **OWNERSHIP OF MATERIALS:**

Materials and information contained in this document are **Simport Scientific**'s copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.

Issue Date: 2017-06-15