

**Product Name:** Dilution vial**Catalogue No.:** V130

This document replaces any previous version

**1. Product Description:**

- **25 ml Dilution vial:** Non sterile, 25 ml vials with caps in trays of 50 units

**2. Packaging:**

- **Case:** 20 trays of 50 vials / 1000 unit per case

**3. Product Specifications:**

- **Material:**
  - Tube: Polystyrene
  - Cap: High Density Polyethylene
  - Tray: Polystyrene
- **Temperature range:** -80°C to +80°C
- **Non-autoclavable**

**4. Standards and Conformity:**

- **ISO 2859-1:** Sampling and inspection procedures

**5. Quality Assurance:**

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

**6. Traceability:**

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

**7. Storage Conditions:**

- Store at room temperature in normal warehouse conditions.

**OWNERSHIP OF MATERIALS:**

Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.