

**Product Name:** Combi-Box™

**Catalogue No.:** T327

This document replaces any previous version

## 1. Product Description:

- **Combi-Box:** Non-sterile, can be used as a storage rack and as a workstation
  
- **Available colors:**
  - White base with transparent cover
  
- **Compatible with:**
  - T327-1 Combi-Rack™
  - Will accept all 96 and 384 well PCR plates

## 2. Packaging:

- **Case:** 5 units per case

## 3. Product Specifications:

- **Material:**
  - PVC
- **Temperature range:** -30°C to +50°C
- **Autoclavable:** No

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions;

## 8. Recommended Use:

Rack for holding tubes

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*