

**Product Name:** DissecTable™ Jr

**Catalogue No.:** M618

## 1. Product Description:

- **M618 Dissecting board:** A smaller dissecting board, perfect for smaller counter area. Printed on both sides. Large and small cutting areas including dimensional scales in inches and centimeters.

## 2. Packaging:

- **Case:** 1 individually wrapped unit per case

## 3. Product Specifications:

- **Material:**
  - High-Density Polyethylene
- **Temperature range:**
  - Base -80° C TO + 121° C
- **Dimensions:** 13 x 11 x 0.5 in.

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging or product

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*