


# R225 Reagent Reservoir Family

*Simport*



**Simport**<sup>®</sup>  
Since 1975 *Scientific inc.*

*A family owned  Canadian company*

[simport.com](http://simport.com)



## R225 Reagent Reservoir for Automation

Packing size : Qty/Pk 10 - Qty/Cs 50

Made from polypropylene.

For labs involved with fraction collecting, compound storage, plate replication, and mixing applications. Suitable for a range of liquid handling configurations from 8 or 12 channel pipettes to 96 or 384 pipette heads. Ideal for small volume robotic tips used in high throughput instruments. Choose from a single reservoir or partitioned models for several liquids.


Conforms to ANSI/SLAS standard dimensions, they are compatible with virtually any automated liquid handler, stacker, and manual procedure.

Made from natural virgin polypropylene, providing a low binding surface to prevent samples from sticking to the sidewalls during elution and is chemically inert for combinatorial chemistry applications.


### Features

- Robot friendly, designed for simple, economical handling of reagents in automated systems
- Bottom designed to enable funneling of liquid, maximizing reagent recovery with pipettes
- Low-profile
- Chemical and heat resistant, can be stored at temperatures as low as -80 °C
- Autoclavable
- The cover cannot go in the autoclave
- Can be gamma irradiated
- Certified DNase/RNase and Pyrogen Free

### Cover





No Cat. **R225-15**  
Material Medical grade polystyrene  
Size mm 127.8 x 85 x 9  
Qty/Pk 10  
Qty/Cs 300





1 Well

Cat. No. **R225-1**  
Volume 190 ml  
Size mm 127.8 x 85.5 x 31.4





8 Wells

Cat. No. **R225-3**  
Total Volume 180 ml  
Per Channel 22 ml  
Size mm 127.8 x 85.5 x 31.4





12 Wells

Cat. No. **R225-5**  
Total Volume 180 ml  
Per Channel 15 ml  
Size mm 127.8 x 85.5 x 31.4



96 Wells

Cat. No. **R225-7**  
Volume 195 ml  
Size mm 127.8 x 85.5 x 31.4



384 Wells

Cat. No. **R225-9**  
Volume 185 ml  
Size mm 127.8 x 85.5 x 31.4

