

Technical Data Sheet

Product Name: Cryovial®

Catalogue No.: T210-X Series

This document replaces any previous version

1. Product Description:

- **Assembled External Thread Cryogenic Vial with Silicone O-Ring Seal and Graduations:**
 - Sterile, Disposable Cryogenic Vial
- **Available Options:**
 - T210-2A - Cryogenic 2D Vials Tubes & Closures (500/Cs)
 - T210-2A81B – Pre-Racked Tubes (10/Cs)
- **Compatible with:**
 - T214-281B Storage Box for the T210 Cryogenic 2D Vials
 - Color Coding Cap Insert T312

2. Packaging:

- **Case:** 5 bags if 100 units / 500 units per case

3. Product Specifications:

- **Material:**
 - Tube: Polypropylene
 - Cap: Polypropylene
 - O-Ring: Silicone
- **Temperature Range:** -196°C to +121°C
- **Autoclavable:** +121°C for up to 30 minutes
- **Leak-proof:** Tested in vacuum at 71.3cm Hg (95 kPa)
- **Liquid Nitrogen (LN²) resistance:** Gas phase of LN² only
- **MEA:** Product has passed MEA testing

4. Standards and Conformity:

- **Radiation & Dosage:** Sterilizer certifies a Gamma SAL of 10⁻³; between 6.5 kGy & 13.5 kGy
- **RN/DNase, DNA & PCR free:** Certified RNase, DNase, DNA & PCR Inhibitor Free
- **Pyrogen-Free:** Certified Acceptance Criteria per USP <85> Bacterial Endotoxins Test met [λ = 0.06 EU/mL; results < 0.06 EU/mL or < 21 EU/device]
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
Silicone O-Ring conforms to FDA 21 CFR 177.2600
- **USP:** Plastic Conforms to USP Class VI
- **CE:** Product is CE Marked
- **CONEG / RoHS:** Suppliers certify resin and O-Ring are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Supplier certifies the resin is in conformity to REACH standards
- **EU Reg. 2019/1021 (POPs):** Supplier certifies the resin is in conformity with EU 2019/1021 (POP)
- **Toxic Substance Control Act:** Supplier certifies the resin is in conformity with US TSCA
- **PFAS:** Supplier certifies the resin as PFAS-Free
- **BSE/TSE:** Supplier certifies the resin as BSE/TSE-Free
- **Bisphenols:** Supplier certifies the resin as free from Bisphenols & their compounds
- **Phthalates:** Supplier certifies the resin as Phthalate Free
- **LATEX:** Supplier certifies the resin as Latex Free
- **Conflict Minerals Reporting:** Supplier certifies the resin is in conformity with CMR standards

Technical Data Sheet

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Volume measurements
- Closure verification
- Leak proof testing in vacuum at 71.3cm Hg
- Gas phase of Liquid Nitrogen resistance

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging which could compromise the product sterility

8. Recommended Use:

- Designed for storing biological material at temperatures as low as -196°C
- Not recommended in the liquid phase of LN². Disposables for single use only.
- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions
- Should be used only in the gas phase of Liquid Nitrogen (LN²)

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*

Disclaimer:

The information contained in this document shall not constitute a representation or warranty of the suitability for any specific purpose. To the best of our knowledge, the information contained is reliable and true based on our internal and external testing as well as statements supplied to us by third party suppliers.